


5F LONG

mc-peripheral

Selfexpandable Nitinol Stent

Open cell design

- Flexible stent design
- Minimal foreshortening



Braided outer sheath



- Minimal force for stent release
- Extremely kink resistant and flexible
- Excellent pushability for cross-over

Application device

- High flexibility
- Radiopaque platinum markers on the inner catheter
- Two different working lengths: 80 cm / 120 cm
- For use with .018 / .021 inch guide wire
- Simple to use

Stent visibility

- Integrated laser-welded tantalum markers at proximal and distal end



medicut
Stent Technology GmbH

Product Specifications

Application Device 5F / 80 cm

Nominal length mm	Stent Ø 5 mm	Stent Ø 6 mm	Stent Ø 7 mm	Stent Ø 8 mm	Stent Ø 9 mm	Stent Ø 10 mm	Order-No.
120	0580-5120	0580-6120	0580-7120	0580-8120	0580-9120	0580-1120	
150	0580-5150	0580-6150	0580-7150	0580-8150			
180	0580-5180	0580-6180	0580-7180				
200	0580-5200	0580-6200	0580-7200				

Application Device 5F / 120 cm

Nominal length mm	Stent Ø 5 mm	Stent Ø 6 mm	Stent Ø 7 mm	Stent Ø 8 mm	Stent Ø 9 mm	Stent Ø 10 mm	Order-No.
120	0513-5120	0513-6120	0513-7120	0513-8120	0513-9120	0513-1120	
150	0513-5150	0513-6150	0513-7150	0513-8150			
180	0513-5180	0513-6180	0513-7180				
200	0513-5200	0513-6200	0513-7200				

Recommended introducer sheath 5F / Maximum guide wire .018 / .021 inch

Distributed by: