


6F LONG

mc-peripheral

Selfexpandable Nitinol Stent

Open cell design

- Flexible stent design
- Minimal foreshortening



Braided outer sheath



- Minimal force for stent release
- Extremely kink resistant and flexible
- Excellent pushability for cross-over

Application device

- High flexibility
- Radiopaque platinum markers on the inner catheter
- Two different working lengths: 80 cm / 130 cm
- For use with .035 inch guide wire
- Simple to use

Stent visibility

- Integrated laser-welded tantalum markers at proximal and distal end



medicut
Stent Technology GmbH

Product Specifications

Application Device 6F / 80 cm

Nominal length mm	Stent Ø 5 mm	Stent Ø 6 mm	Stent Ø 7 mm	Stent Ø 8 mm	Stent Ø 9 mm	Stent Ø 10 mm	Order-No.
120	0680-5120	0680-6120	0680-7120	0680-8120	0680-9120	0680-1120	
150	0680-5150	0680-6150	0680-7150	0680-8150			
180	0680-5180	0680-6180	0680-7180				
200	0680-5200	0680-6200	0680-7200				

Application Device 6F / 130 cm

Nominal length mm	Stent Ø 5 mm	Stent Ø 6 mm	Stent Ø 7 mm	Stent Ø 8 mm	Stent Ø 9 mm	Stent Ø 10 mm	Order-No.
120	0613-5120	0613-6120	0613-7120	0613-8120	0613-9120	0613-1120	
150	0613-5150	0613-6150	0613-7150	0613-8150			
180	0613-5180	0613-6180	0613-7180				
200	0613-5200	0613-6200	0613-7200				

Recommended introducer sheath 6F / Maximum guide wire .035 inch

Distributed by: _____
