

**.014 inch**

# mc-PTA

Balloon Catheter



Semi-compliant balloon, diameters from 2 to 4 mm and lengths from 40 to 150mm


# medicut

Stent Technology GmbH

# mc-PTA Balloon Catheter

**.014 inch**

- Coaxial shaft design for short inflation and deflation time
- Flexible kink-resistant shaft
- Short conical tip
- Semi-compliant balloon
- Markers define the length from outer edge to outer edge
- Hydrophilic coated
- Low entry balloon design



Order Code usable shaft length		Ballon Ø Ø / mm	Ballon length cm	Shaft Ø F	Recommended introducer size F	Nominal pressure bar	Rated burst pressure bar
120cm	150cm						
1412-2004	1415-2004	2	4	3,9	4	6	16
1412-2008	1415-2008	2	8	3,9	4	6	16
1412-2012	1415-2012	2	12	3,9	4	6	16
1412-2015	1415-2015	2	15	3,9	4	6	16
1412-2504	1415-2504	2,5	4	3,9	4	6	16
1412-2508	1415-2508	2,5	8	3,9	4	6	16
1412-2512	1415-2512	2,5	12	3,9	4	6	16
1412-2515	1415-2515	2,5	15	3,9	4	6	16
1412-3004	1415-3004	3	4	3,9	4	6	14
1412-3008	1415-3008	3	8	3,9	4	6	14
1412-3012	1415-3012	3	12	3,9	4	6	14
1412-3015	1415-3015	3	15	3,9	4	6	14
1412-3504	1415-3504	3,5	4	3,9	4	6	14
1412-3508	1415-3508	3,5	8	3,9	4	6	14
1412-3512	1415-3512	3,5	12	3,9	4	6	14
1412-3515	1415-3515	3,5	15	3,9	4	6	14
1412-4004	1415-4004	4	4	3,9	4	6	14
1412-4008	1415-4008	4	8	3,9	4	6	14
1412-4012	1415-4012	4	12	3,9	4	6	14
1412-4015	1415-4015	4	15	3,9	4	6	14

Adapted to .014 inch guide wire

Distributed by:

**medicut** Stent Technology GmbH

Wilhelm-Becker-Straße 11a · D-75179 Pforzheim / Germany

Fon +49 (0) 7231-280 100 · Fax +49 (0) 7231-280 1025 · info@medi-cut.de · www.medi-cut.de