



.035 inch

mc-PTA

Balloon Catheter



Semi-compliant balloon, diameters from 4 to 8 mm and lengths from 20 to 150mm


medicut

Stent Technology GmbH

mc-PTA Balloon Catheter

.035 inch


- Bilumen design
- Flexible kink-resistant shaft
- Short conical tip
- Semi-compliant balloon
- Markers define the length from outer edge to outer edge
- Hydrophilic coated



Order Code usable shaft length		Ballon Ø Ø / mm	Ballon length cm	Shaft Ø F	Recommended introducer size F	Nominal pressure bar	Rated burst pressure bar
80cm	135cm						
3508-4002	3513-4002	4	2	5	5	6	16
3508-4004	3513-4004	4	4	5	5	6	16
3508-4006	3513-4006	4	6	5	5	6	16
3508-4008	3513-4008	4	8	5	5	6	16
3508-4010	3513-4010	4	10	5	5	6	16
3508-4012	3513-4012	4	12	5	6	6	16
3508-4015	3513-4015	4	15	5	6	6	16
3508-5002	3513-5002	5	2	5	5	6	16
3508-5004	3513-5004	5	4	5	5	6	16
3508-5006	3513-5006	5	6	5	5	6	16
3508-5008	3513-5008	5	8	5	5	6	14
3508-5010	3513-5010	5	10	5	5	6	14
3508-5012	3513-5012	5	12	5	5	6	14
3508-5015	3513-5015	5	15	5	5	6	14
3508-6002	3513-6002	6	2	5	5	6	16
3508-6004	3513-6004	6	4	5	5	6	16
3508-6006	3513-6006	6	6	5	5	6	16
3508-6008	3513-6008	6	8	5	5	6	14
3508-6010	3513-6010	6	10	5	5	6	14
3508-6012	3513-6012	6	12	5	6	6	14
3508-6015	3513-6015	6	15	5	6	6	14
3508-7002	3513-7002	7	2	5	6	6	14
3508-7004	3513-7004	7	4	5	6	6	14
3508-7006	3513-7006	7	6	5	6	6	14
3508-7008	3513-7008	7	8	5	6	6	12
3508-7010	3513-7010	7	10	5	6	6	12
3508-7012	3513-7012	7	12	5	6	6	12
3508-7015	3513-7015	7	15	5	6	6	12
3508-8002	3513-8002	8	2	5	6	6	14
3508-8004	3513-8004	8	4	5	6	6	14
3508-8006	3513-8006	8	6	5	6	6	14
3508-8008	3513-8008	8	8	5	6	6	12
3508-8010	3513-8010	8	10	5	6	6	12
3508-8012	3513-8012	8	12	5	6	6	12
3508-8015	3513-8015	8	15	5	6	6	12

Adapted to .035 inch guide wire

Distributed by:



medicut Stent Technology GmbH

Wilhelm-Becker-Straße 11a · D-75179 Pforzheim / Germany

Fon +49 (0) 7231-280 100 · Fax +49 (0) 7231-280 1025 · info@medi-cut.de · www.medi-cut.de