

V+Pad™ Non-Invasive Hemostasis Pad

The V+Pad™ Non-Invasive Hemostasis Pad is a hydrophilic wound dressing intended for use in the local management of bleeding wounds such as vascular access sites, percutaneous catheters or tubes, surgical debridement and lacerations. V+Pad™ promotes rapid control of bleeding in patients on anticoagulation therapy.

- Controls persistent bleeding at central venous and arterial catheter sites, and may be cut for easy placement around vascular access catheters
- Double-stitched, woven gauze composition
- D-glucosamine enriched gauze fibers
- Use in conjunction to manual compression



V+Pad™ Non-Invasive Hemostasis Pad (Box of 10)

CATALOG NO.	DESCRIPTION
3935	V+Pad™ Non-Invasive Hemostasis Pad, approx. 4 cm x 4 cm