

# Fine Needle Aspiration (FNA) Biopsy Needles

Our traditional biopsy needles harvest aspirant for cytological evaluations or specimens for histological evaluations.

## Disposable Soft Tissue Biopsy Needles:

- Each needle has numerically ordered centimeter markings on the outer cannula to facilitate precise depth placement
- Needles feature an adjustable needle stop which allows the user to restrict forward movement, localizing the needle tip to the biopsy site
- Color coded stylet hubs indicate gauge size of needles
- Available in a variety of gauge sizes and centimeter lengths
- Available with an echogenic tip on the outer cannula to promote accurate placement under ultrasound guidance

### Chiba Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MCN1803	18G x 9cm
MCN1803US	18G x 9cm, echogenic tip
MCN1806	18G x 15 cm
MCN1806US	18G x 15 cm, echogenic tip
MCN1808	18G x 20 cm
MCN1808US	18G x 20 cm, echogenic tip
MCN2003	20G x 9 cm
MCN2003US	20G x 9 cm, echogenic tip
MCN2006	20G x 15 cm
MCN2006US	20G x 15 cm, echogenic tip
MCN2008	20G x 20 cm
MCN2008US	20G x 20 cm, echogenic tip
MCN2203	22G x 9 cm
MCN2203US	22G x 9 cm, echogenic tip
MCN2206	22G x 15 cm
MCN2206US	22G x 15 cm, echogenic tip
MCN2208	22G x 20 cm
MCN2208US	22G x 20 cm, echogenic tip



Chiba

### Spinal Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MSPN1803	18G x 9 cm
MSPN1806	18G x 15 cm
MSPN2003	20G x 9 cm
MSPN2003US	20G x 9 cm, echogenic tip
MSPN2006	20G x 15 cm
MSPN2006US	20G x 15 cm, echogenic tip
MSPN2203	22G x 9 cm
MSPN2203US	22G x 9 cm, echogenic tip
MSPN2206	22G x 15 cm
MSPN2206US	22G x 15 cm, echogenic tip
MSPN2208	22G x 20 cm



Spinal

For more information, please call the Angiotech Customer Care Center at 877 991-1110 (USA and Canada only) or 732 626-6466, option 1  
 • Email [interventional@angio.com](mailto:interventional@angio.com) • Visit [www.AngioEduPRO.com](http://www.AngioEduPRO.com), our educational website for healthcare professionals

# Fine Needle Aspiration (FNA) Biopsy Needles

## Franseen Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MFN1803	18G x 9 cm
MFN1803US	18G x 9 cm, echogenic tip
MFN1806	18G x 15 cm
MFN1806US	18G x 15 cm, echogenic tip
MFN1808	18G x 20 cm
MFN1808US	18G x 20 cm, echogenic tip
MFN2003	20G x 9 cm
MFN2003US	20G x 9 cm, echogenic tip
MFN2006	20G x 15 cm
MFN2006US	20G x 15 cm, echogenic tip
MFN2008	20G x 20 cm
MFN2008US	20G x 20 cm, echogenic tip
MFN2203	22G x 9 cm
MFN2203US	22G x 9 cm, echogenic tip
MFN2206	22G x 15 cm
MFN2206US	22G x 15 cm, echogenic tip
MFN2208	22G x 20 cm
MFN2208US	22G x 20 cm, echogenic tip

## Greene Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MGN2206X	22G x 15 cm incl. 19G x 10 cm introducer, echogenic tip

## Westcott Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MWN2003	20G x 9 cm
MWN2003US	20G x 9 cm, echogenic tip
MWN2006	20G x 15 cm
MWN2006US	20G x 15 cm, echogenic tip
MWN2008	20G x 20 cm
MWN2008US	20G x 20 cm, echogenic tip
MWN2203	22G x 9 cm
MWN2203US	22G x 9 cm, echogenic tip
MWN2206	22G x 15 cm
MWN2206US	22G x 15 cm, echogenic tip
MWN2208	22G x 20 cm
MWN2208US	22G x 20 cm, echogenic tip

## Locking Syringe (Box of 10) (not illustrated)

CATALOG NO.	DESCRIPTION
MLS1001	10 ml



Franseen



Greene



Westcott