

Osty-Core™ Bone Biopsy Needle

The unique tip geometry of the Osty-Core™ Bone Biopsy Needle is designed to provide exceptional needle entry and coring characteristics without distorting the biopsy specimen.

- The ergonomic handle features a twist-lock design that stabilizes the needle during placement
- A luer lock connector on the handle for secure syringe attachment
- An adjustable, removable needle stop safely controls depth of penetration
- Can be used with co-axial introducer needles



OBN

Osty-Core™ Bone Biopsy Needle (Box of 10)

CATALOG NO.	DESCRIPTION
OBN16075	16G x 7.5 cm
OBN16125	16G x 12.5 cm, optional co-axial needle MCXS1610AX
OBN16165	16G x 16.5 cm, optional co-axial needle MCXS1616AX

Optional Co-axial Introducer Needles (sold separately) (Box of 10)

CATALOG NO.	DESCRIPTION
MCXS1610AX	15G x 4.6 cm, co-axial to OBN16125
MCXS1616AX	15G x 10.6 cm, co-axial to OBN16165