

Bone Marrow Aspiration Needles

Safe aspirations at the sternum can be performed by adjusting the needle stop to control the depth of penetration. The needle stop can be removed if aspirating at the iliac crest.

Features include:

- Luer lock connector on the handle for secure syringe attachment
- An adjustable, removable needle stop safely controls depth of penetration

Bone Marrow Aspiration Needle (Box of 10)

CATALOG NO.	DESCRIPTION
DBMNI1501	15G x 2.688" max.
DBMNI1504	15G x 4" max.
DBMNI1601	16G x 2.688" max.



Adjustable needle stop

For more information, please call the Angiotech Customer Care Center at 877 991-1110 (USA and Canada only) or 732 626-6466, option 1
• Email interventional@angio.com • Visit www.AngioEduPRO.com, our educational website for healthcare professionals