

Lunderquist Stainless Steel Guidewires

The Lunderquist Guidewire is an extremely stiff guidewire used to provide extra support during non-vascular procedures. Flexible tip for access to delicate or difficult areas.

- Lunderquist guidewire for firm support.
- Available in a variety of lengths and tip configurations.

Lunderquist Stainless Steel Guidewire • Straight **Box of 10**

1051 28 082 .028" (0.71 mm) x 82 cm, 7.5 cm flexible coil

1051 35 082 .035" (0.89 mm) x 82 cm, 7.5 cm flexible coil

1051 35 120 .035" (0.89 mm) x 120 cm, 7.5 cm flexible coil

Lunderquist Stainless Steel Guidewire • 3 mm J Curve **Box of 10**

1052 28 082 .028" (0.71 mm) x 82 cm, 7.5 cm flexible coil

1052 35 082 .035" (0.89 mm) x 82 cm, 7.5 cm flexible coil

1052 35 120 .035" (0.89 mm) x 120 cm, 7.5 cm flexible coil

