

POINTER® Nitinol Guidewire

Flexibility and stability, strength and malleability all in one wire

– available in two sizes: .018" and .020"

- Excellent visibility through highly radiopaque coil tip.
- Smooth transition from nitinol core to hydrophilic coated radiopaque tip provides excellent manoeuvrability.
- HYDROFEEL® hydrophilic coating on tip for less traumatic path finding/advancement.
- PTFE coated nitinol shaft for minimal friction during introduction and manipulation.
- Nitinol shaft construction provides 1:1 torque control and kink-resistance for negotiating through tortuous anatomy.

POINTER® Nitinol Guidewire • 4 cm flexible coil *Box of 5*


1804 18 210 .018" (0.46 mm) x 210 cm

1804 18 260 .018" (0.46 mm) x 260 cm

1804 18 300 .018" (0.46 mm) x 300 cm

1804 20 210 .020" (0.49 mm) x 210 cm

1804 20 300 .020" (0.49 mm) x 300 cm



POINTER® Nitinol Guidewire • 6 cm flexible coil *Box of 5*

1806 20 210 .020" (0.49 mm) x 210 cm

1806 20 300 .020" (0.49 mm) x 300 cm