

Access Guidewires

Access Guidewires

- High torque-control.
- Smooth transition from rigid core to radiopaque coil tip provides excellent manoeuvrability.
- Nitinol offers kink resistance.
- Stainless steel for maximum back-up.

.018 Nitinol Access Guidewire

Box of 5

1726 18 080 .018" (0.46 mm) x 80 cm, 6.5 cm radiopaque coil tip

.018 Access Guidewire • Stainless Steel Shaft

Box of 5

1070 18 080 .018" (0.46 mm) x 80 cm, 6.5 cm radiopaque coil tip

