

Screw Dilators with Introducer Split Sheath


A single screw dilator replaces several conventional dilators. The screw tip penetrates tissue more thoroughly, enabling dilation up to 14F using only one dilator. The split sheath offers free passage for the catheter and is simple to remove by peeling it apart.

- Blunt dissection.
- Introducer split sheath for protection of the percutaneous tract.
- Tissue dilation up to 14F using only one dilator.
- Radiopaque for fluoroscopic visualization screw dilators with introducer split sheaths.



Screw Dilators with Introducer Split Sheath **Box of 10**

- 6407 38 025** 7F x 25 cm, accepts guidewire .038" (0.97 mm)
- 6408 38 025** 8F x 25 cm, accepts guidewire .038" (0.97 mm)
- 6409 38 025** 9F x 25 cm, accepts guidewire .038" (0.97 mm)
- 6410 38 025** 10F x 25 cm, accepts guidewire .038" (0.97 mm)
- 6412 38 025** 12F x 25 cm, accepts guidewire .038" (0.97 mm)
- 6414 38 025** 14F x 25 cm, accepts guidewire .038" (0.97 mm)
- 6416 47 025** 16F x 25 cm, accepts guidewire .047" (1.19 mm)



64xx