

# SKATER® Introducer Sets

The SKATER® Introducer Set provides easy, accurate and atraumatic placement of .035" or .038" guidewires in non-vascular procedures.

- Smooth transition between the 4 F dilator and the 6 F sheath.
- Embedded gold ring marker enhances visualization for accurate fluoroscopic placement.
- 21G needle minimizes tissue trauma.
- .018" guidewire with malleable coil tip for easy advancement.



6545

## SKATER® Introducer Set

For introduction of .035" or .038" guidewire **Box of 5**

- 6515 06 000** Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm stainless steel access guidewire, 6F x 18 cm co-axial dilator
- 6515 06 300** Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access guidewire, 6F x 18 cm co-axial dilator
- 6555 06 000** Set, 21G (0.8 mm) x 15 cm chiba needle, .018" (0.46 mm) x 80 cm stainless steel access guidewire, 6F x 18 cm co-axial dilator
- 6555 06 300** Set, 21G (0.8 mm) x 15 cm chiba needle, .018" (0.46 mm) x 80 cm nitinol access guidewire, 6F x 18 cm co-axial dilator

## SKATER® Introducer Set • Biliary Access

For introduction of .035" or .038" guidewire **Box of 5**

- 6514 06 000** Biliary Access Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm stainless steel access guidewire, 6F x 26 cm co-axial dilator
- 6514 06 300** Biliary Access Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access guidewire, 6F x 26 cm co-axial dilator
- 6554 06 000** Biliary Access Set, 21G (0.8 mm) x 15 cm chiba needle, .018" (0.46 mm) x 80 cm stainless steel access guidewire, 6F x 26 cm co-axial dilator
- 6554 06 300** Biliary Access Set, 21G (0.8 mm) x 15 cm chiba needle, .018" (0.46 mm) x 80 cm nitinol access guidewire, 6F x 26 cm co-axial dilator

## SKATER® Introducer Set • Complete

**Box of 5**

- 6545 06 200** Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access guidewire, 6F x 18 cm co-axial dilator, .035" (0.89 mm) x 150 cm Amplatz guidewire

## SKATER® Introducer Biliary Set • Complete

**Box of 5**

- 6543 06 200** Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access guidewire, 6F x 26 cm co-axial dilator, .035" (0.89 mm) x 150 cm Amplatz guidewire

Drainage



Smooth transition between .018" guide wire and SKATER® Introducer



The white sheath accepts both a .018" and a .038" guidewire