

SKATER™ Centesis Catheters

Designed for percutaneous fluid aspirations and small volume drainage procedures.

- Vent plug allows air to exit catheter during insertion and protects clinician from fluid leakage
- Luer locked introducer needle for secure one-handed placement and hollow stylet needle design permits aspiration of fluid prior to stylet removal
- Tapered catheter tip ensures good transition between catheter and needle to facilitate smooth insertion
- Four side holes on the distal tip of the catheter facilitate drainage within small cavities



7003

SKATER™ Centesis Catheters (Box of 10)

CATALOG NO.	DESCRIPTION
700304007	4F x 7 cm
700304010	4F x 10 cm
700304015	4F x 15 cm
700305007	5F x 7 cm
700305010	5F x 10 cm
700305015	5F x 15 cm
700305020	5F x 20 cm