

Drainage Bags

350 ml sterile drainage bag with flocked back and bottom outlet. Should be used with a male to male connecting tube e.g. 4040 40 xxx for secure connection to the drainage catheter.

- Flocked back reduces sweating.
- Minimal patient discomfort.
- Anti-reflux valve to prevent retrograde contamination.
- Acceptable size – can be hidden beneath clothes.
- High patient acceptance.

8500



8500 07 350 350 ml drainage bag – female.

Box of 50

600 ml sterile drainage bag with bottom outlet and adjustable drain tube length.

- Soft kink-resistant tube can be cut and fit to the appropriate length.
- Twist type valve for quick and easy emptying of drainage bag.
- Anti-reflux valve to prevent retrograde contamination.
- Fully adjustable soft belt.

DBAG



DBAG600 600 ml drainage bag with adjustable drain tube – male luer lock.

Box of 20