

# Galactography Kits

Galactography kits are designed to perform safe and effective contrast media enhanced radiography of the mammary ducts.


- Dilators and injection cannula have blunt ends to minimize puncture of the duct walls
- Dilators are color coded for easy diameter identification
- Injection cannula is available with a curved tip to facilitate insertion

## Galactography Kit (Box of 10)

CATALOG NO.	DESCRIPTION
DGK002SX	30G w/straight cannula
DGK002X	30G w/curved cannula
RTD001X	24G w/curved cannula



DGK002X



RTD001X

For more information, please call the Angiotech Customer Care Center at 877 991-1110 (USA and Canada only) or 732 626-6466, option 1  
 • Email [interventional@angiotech.com](mailto:interventional@angiotech.com) • Visit [www.AngioEduPRO.com](http://www.AngioEduPRO.com), our educational website for healthcare professionals