

mediform

SPECIÁLNÍ  
ZDRAVOTNICKÝ  
MATERIÁL

Joline®

# Allevo

Kyphoplasty System



**Balloon in perfection**

Catheter specialists for more than 10 years

Made in  
Germany

# Allevo

## Kyphoplasty System

### Allevo Kyphoplasty System - Advanced stabilization & pain relief for spinal fractures

Joline's minimally-invasive balloon catheters provide a complete system for treating vertebral fractures. Our unique balloon design allows our balloon catheters to be precisely positioned in the vertebra, thereby creating the best conditions for stabilizing the spine. Years of experience and precision engineering combine to provide numerous advantages and guarantee users the highest degree of safety and reliability.


Joline — Close to the customer.

Basic  
Kit

Smart  
Kit

Best  
Kit

The convenience of having the most  
suitable Kit for every Kyphoplasty treatment



Color-coding makes our  
instruments quick & easy to identify



S5001	<b>One Pedicle Kit</b>	Balloon Ø 16 mm; Balloon Length = 10mm	<b>Kit contains:</b> 1 x Vertebra Access Device 1x each Guide Pin, blunt, sharp 1 x Vertebra Introducer Device 1 x Vertebra Biopsy Device 1 x Bone Drill 1 x appropriate balloon catheter 1 x Vacuum Syringe 1 x Inflation Device 3 x Vertebra Filler Device 1 x Bone cement
S5002	<b>One Pedicle Kit</b>	Balloon Ø 16 mm; Balloon Length = 16mm	
S5003	<b>One Pedicle Kit</b>	Balloon Ø 16 mm; Balloon Length = 22mm	

## Allevo

Order information

**Individual components**  
The Allevo kyphoplasty system comprises an extensive selection of individual instruments



**Kits** Joline offers different system combinations with all kits customized to the exact need of the user



On the right  
side you can find  
the kits.

mediform

SPECIÁLNÍ  
ZDRAVOTNICKÝ  
MATERIÁL

MEDIFORM, spol. s r.o., Olšová 1, CZ Brno 637 00  
TELEFON: +420 541 221 392 / FAX: +420 541 220 736  
E-MAIL: objednavky@mediform.cz / www.mediform.cz