


# Sclerotherapy Injection Needles

**Used for endoscopic injection into esophageal varices above esophageal junction**



Sharp-Bevel-Needle with JET-Pressure effect metal protection

Especially developed needle to provide secure firm grip with inner tube and prevent possible leakage from the joint of inner tube and needle.

All tubes are made of PTFE (Teflon)

Art. No.	Product name	Needle gauge	Needle length	Working length	Working Channel Ø
SGN-2321	Sclerotherapy Needle	21G	7 mm	170 cm	2.8 mm
SGN-2321-C	Sclerotherapy Needle	21G	7 mm	230 cm	2.8 mm
SGN-2323	Sclerotherapy Needle	23G	7 mm	170 cm	2.8 mm
SGN-2323-C	Sclerotherapy Needle	23G	7 mm	230 cm	2.8 mm
SGN-2324	Sclerotherapy Needle	24G	4 mm	170 cm	2.8 mm
SGN-2324-C	Sclerotherapy Needle	24G	4 mm	230 cm	2.8 mm

C = Colono