


# Loop Basket

**Used for Endoscopic Retrieval of foreign bodies from the Gastrointestinal Tract and to retrieve polyp through colono scope**

**Shape: oval with flush port**

**Sheath: PTFE**



Flat Basket available on request

Art. No.	Product name	Loop Size	Working Length	Working Channel Ø
SFBH-LB	Loop Basket	60 mm	170 cm	2.8 mm
SFBH-LB-C	Loop Basket	60 mm	230 cm	2.8 mm
SFBH-LB-40	Loop Basket	40 mm	170 cm	2.8 mm
SFBH-LB-C-40	Loop Basket	40 mm	230 cm	2.8 mm
SFBH-LB-90	Loop Basket	90 mm	170 cm	2.8 mm
SFBH-LB-C-90	Loop Basket	90 mm	230 cm	2.8 mm