


# Biliary Stone Extraction Balloon Catheter

Wire guided, triple lumen

Used for endoscopic extraction  
of Biliary Stones.



Art. No.	Catheter Diameter mm/FR	Functional Catheter Length	Diameter of Inflated Balloon	Balloon Max. Gas Inflation	Balloon Max. Fluid Injection	Balloon Length	Recommended Guidewire	Syringe
S-SEBC-7	2.35/7	200 cm	14-16 mm	4.5 cm <sup>2</sup>	3 ml	16 mm	0.035"	5 cc