



Ureteral Stents

- **Standard Double-J**
- **Multiloop Stents**



New pusher with x-ray marking



Ureteral Stents

Standard Double-J and Multiloop Stents

1 set consists of Double-J stent with suture and pusher with x-ray marking, with or without guidewire, Stent made of aliphatic polyurethane (Tecoflex®) Guidewire stainless steel, PTFE coated, one side stiff, one side flexible
150 cm, 0.032", 0.035" or 0.038"

Closed tip

Open tip

Art. No.	Art.No.	Size (Fr/CH)	Length (cm)	Guidewire
SS-4716	SOT-4716	4.7	16	0.032"
SS-4724	SOT-4724	4.7	24	0.032"
SS-4726	SOT-4726	4.7	26	0.032"
SS-4728	SOT-4728	4.7	28	0.032"
SS-4730	SOT-4730	4.7	30	0.032"
SS-47 ML+ML	SOT-47 ML+ML	4.7	24 -32	0.032"
SS-516	SOT-516	5	16	0.035"
SS-524	SOT-524	5	24	0.035"
SS-526	SOT-526	5	26	0.035"
SS-528	SOT-528	5	28	0.035"
SS-530	SOT-530	5	30	0.035"
SS-5 ML+ML	SOT-5 ML+ML	5	2x Multiloop	0.035"
SS-616	SOT-616	6	16	0.035"
SS-624	SOT-624	6	24	0.035"
SS-626	SOT-626	6	26	0.035"
SS-628	SOT-628	6	28	0.035"
SS-630	SOT-630	6	30	0.035"
SS-6 ML+ML	SOT-6 ML+ML	6	2x Multiloop	0.035"
SS-716	SOT-716	7	16	0.038"
SS-724	SOT-724	7	24	0.038"
SS-726	SOT-726	7	26	0.038"
SS-728	SOT-728	7	28	0.038"
SS-730	SOT-730	7	30	0.038"
SS-7 ML+ML	SOT-7 ML+ML	7	2x Multiloop	0.038"
SS-826	SOT-826	8	26	0.038"
SS-828	SOT-828	8	28	0.038"
SS-830	SOT-830	8	30	0.038"
SS-8 ML+ML	SOT-8 ML+ML	8	2x Multiloop	0.038"

other lengths available

If you like the set with guidewire please add +G (Example SS-4716+G)