


Ureteral Stents

- Mono-J Stents
- Blue Stents via Ureterorenoscope



Mono-J Luer-Lock connector

New pusher with x-ray marking

Ureteral Stents

Mono-J Stent

Closed tip	Open tip	Size (Fr/CH)	Length (cm)	Guidewire
MJC-4.770	MJO-4.770	4.7	70	0.032"
MJC-670	MJo-670	6	70	0.035"
MJC-770	MJO-770	7	70	0.038"
MJC-870	MJO-870	8	70	0.038"

"Blue Stents" via Ureterorenoscope with Pusher 70 cm

Closed tip	Open tip	Size (Fr/CH)	Length (cm)	Guidewire
SSU-4716	SOU-4716	4.7	16	0.032"
SSU-4724	SOU-4724	4.7	24	0.032"
SSU-4726	SOU-4726	4.7	26	0.032"
SSU-4728	SOU-4728	4.7	28	0.032"
SSU-4730	SOU-4730	4.7	30	0.032"
SSU-47ML+ML	SOU-47ML+ML	4.7	24-32	0.032"

If you like the set with guidewire please add +G