



Ureteral Stents

Tumor Stent


Reinforced with metal spring
Open tip



Stent reinforced with metal spring



New pusher with x-ray marking



Art. No.	Size (Fr/CH)	Length (cm)	Guidewire
TS 6/24	6	24	0.035"
TS 6/26	6	26	0.035"
TS 6/28	6	28	0.035"
TS 7/24	7	24	0.035"
TS 7/26	7	26	0.035"
TS 7/28	7	28	0.035"

If you like the set with guidewire please add +G (Example TS 6/24+G)