


Ureteral Stents


Anti-Bacteria

Sulfadiazine coated,
a sulfonamide antibiotic

1 set consists of Double J stent (open tip)
and pusher, with (+G) or without guidewire



New pusher with x-ray marking



Art. No.	Size (Fr/CH)	Length (cm)	Guidewire
SOT-4826-AB	4.8	26	0.032"
SOT-4828-AB	4.8	28	0.032"
SOT-626-AB	6	26	0.035"
SOT-628-AB	6	28	0.035"
SOT-726-AB	7	26	0.038"
SOT-728-AB	7	28	0.038"

other sizes available

If you like the set with guidewire please add +G
(Example SOT-4826-AB+G)