

Percutaneous Nephrostomy Drainage Set with 2-way INTEGRAL Balloon Catheter



Percutaneous Nephrostomy Drainage Set with 2-way INTEGRAL Balloon Catheter

A convenient complete set for percutaneous placement of a 2-way Silicone integral balloon catheter. Excellent biocompatibility and drainage.

Nephrostomy Drainage Set

Art. No.	Catheter Size (Fr/CH)	Balloon (cc/ml)	Facial Dilators (Fr/CH)
CH-PCNL 10	10	1.5	6 – 8 – 10
CH-PCNL 12	12	2	8 – 10 – 12
CH-PCNL 14	14	3	10 – 12 – 14
CH-PCNL 16	16	5	10 – 12 – 14 – 16

Nephrostomy Exchange Set

Same as drainage set without puncture needle

CH-PCNL 10-E	10	1.5	6 – 8 – 10
CH-PCNL 12-E	12	2	8 – 10 – 12
CH-PCNL 14-E	14	3	10 – 12 – 14
CH-PCNL 16-E	16	5	10 – 12 – 14 – 16

Nephrostomy Integral Balloon Catheter only

CH-PCNL 10-S	16	1.5
CH-PCNL 12-S	12	2
CH-PCNL 14-S	14	3
CH-PCNL 16-S	16	5
CH-PCNL 18-S	18	5
CH-PCNL 20-S	20	5