


TUR Electrodes


for Prostate Surgery




Sphere Electrode (5 mm Ball)




Cutting Loop for 30° Optic




Necklace Electrode




Cold Knife Electrode



Roller Electrode



Kriloktrode Electrode
(Simultaneous Cutting & Vaporization)



Collins Knife



Spike Roller Electrode

TUR Electrodes, BIPOLAR

For Prostate Surgery

STORZ, 2 Stem, 30° Degrees, 24/26 Fr

Cutting Loop **UBSCP**
Cutting Loop **UBSCP40** (0.40mm wire)

Cutting Loop **UBSC**
Cutting Loop **UBSC30** (0.30mm wire)
Cutting Loop **UBSC40** (0.40mm wire)

Cutting Loop **UBSCDV** (small)
Cutting Loop **UBSCD** (small)
Cutting Loop **UBSCB** (straight)

STORZ, 2 Stem, 30° Degrees, 24/26 Fr

Coagulation Electrode **UBSCL** (pointed)
Coagulation Electrode **UBSCLV** (pointed)

Half moon Vaporization
Electrode, ball-shaped **UBSB**

OLYMPUS, 12° Degrees, 24 Fr

Cutting Loop **UBOC12-S** (small)
Cutting Loop **UBOC12-M** (medium)

OLYMPUS, 30° Degrees, 24 Fr

Cutting Loop **UBOC30-S** (small)
Cutting Loop **UBOC30-M** (medium)

OLYMPUS, 12° and 30° Degrees, 24 Fr

Angeled Loop **UBOCA-S** (small)
Angeled Loop **UBOCA-M** (medium)

OLYMPUS, 12° and 30° Degrees

Roller Electrode **UBOR**
Loop Electrode **UBOC12-D**, 12° large
Loop Electrode **UBOC30-D**, 30° large
Band Electrode **UBOC12-T**, 12° large
Band Electrode **UBOC30-T**, 30° large

OLYMPUS, 12°, 30° and 45° Degrees

Needle,
Angeled Loop
UBOL

Button for Plasma **UBOB**