

# Connecting Tubes – High Pressure

Non-leaking connection to standard luer lock connections

- Standard luer lock connections with or without stopcock.
- Rotational connecting part.

## Connecting Tubes – High Pressure

Box of 25

- 4030 07 050** 50 cm female – male luer lock, max pressure 710 psi
- 4037 11 025** 25 cm female-rotating male luer lock, co-extruded tubing, max pressure 1000 psi
- 4037 11 060** 60 cm female-rotating male luer lock, co-extruded tubing, max pressure 1000 psi
- 4037 11 100** 100 cm female-rotating male luer lock, co-extruded tubing, max pressure 1000 psi
- 4037 11 120** 120 cm female-rotating male luer lock, co-extruded tubing, max pressure 1000 psi
- 8010 00 102** One way stopcock, max pressure 1000 psi
- 8015 00 103** Three-way stopcock/two females/one male luer lock, max pressure 1000 psi




4030 07



4037 11



8010



8015