

# V•Stick™ Vascular Access Set

V•Stick™ Vascular Access Set provides placement of .035" and .038" guidewires into the vascular system using small needle access to reduce puncture site and vessel trauma.

- 21G entry needle with/without echogenic tip
- Radiopaque tip provides excellent manoeuvrability
- V•Stick™ features a smooth transition from inner dilator to the outer sheath
- Available with either .018" stainless steel or nitinol guidewire with platinum tip
- Each set includes a 21G x 7 cm needle, with/without echogenic tip, a .018" x 40 cm stainless steel or nitinol guidewire with platinum tip and a 5F or 4F x 10 cm co-axial dilator/sheath



V•Stick™

## V•Stick™ Vascular Access Set (Box of 10)

CATALOG NO.	DESCRIPTION
391504500E	4F x 10 cm co-axial dilator, 21G x 7 cm needle, .018" x 40 cm stainless steel guidewire with platinum tip
391505500E	5F x 10 cm co-axial dilator, 21G x 7 cm needle, .018" x 40 cm stainless steel guidewire with platinum tip
391504000E	4F x 10 cm co-axial dilator, 21G x 7 cm needle with echogenic tip, .018" x 40 cm stainless steel guidewire with platinum tip
391505000E	5F x 10 cm co-axial dilator, 21G x 7 cm needle with echogenic tip, .018" x 40 cm stainless steel guidewire with platinum tip
391504400E	4F x 10 cm co-axial dilator, 21G x 7 cm needle, .018" x 40 cm nitinol guidewire with platinum tip
391505400E	5F x 10 cm co-axial dilator, 21G x 7 cm needle, .018" x 40 cm nitinol guidewire with platinum tip
391504300E	4F x 10 cm co-axial dilator, 21G x 7 cm needle with echogenic tip, .018" x 40 cm nitinol guidewire with platinum tip
391505300E	5F x 10 cm co-axial dilator, 21G x 7 cm needle with echogenic tip, .018" x 40 cm nitinol guidewire with platinum tip