

MaxiCELL™ Biopsy Needles

MaxiCELL™ Biopsy Needle harvests more aspiration for cytological evaluation than traditional FNA needles.

- Harvests on both forward and backward thrusts of the needle, achieving a greater quantity of harvested cells
- A 30 degree matched ground needle tip geometry helps harvest a cluster of intact cells
- Each needle has numerically ordered centimeter markings on the outer cannula to facilitate depth placement
- Available with an echogenic tip on the outer cannula to promote accurate placement under ultrasound guidance
- Each needle has a unique, adjustable needle stop, which allows you to restrict forward movement, effectively localizing the needle tip to the biopsy site
- Available with co-axial introducer needle



MaxiCELL™

MaxiCELL™ Biopsy Needle (Box of 10)

CATALOG NO.	DESCRIPTION
MMC2202	22G x 5 cm
MMC2203	22G x 9 cm
MMC2206	22G x 15 cm
MMC2202US	22G x 5 cm, echogenic tip
MMC2203US	22G x 9 cm, echogenic tip
MMC2206US	22G x 15 cm, echogenic tip
MMCI2203US	22G x 9 cm, w/ 19G x 5 cm, w/co-axial intr. and echogenic tip
MMCI2206US	22G x 15 cm, w/ 19G x 11.4 cm, w/co-axial intr. and echogenic tip

For more information, please call the Angiotech Customer Care Center at 877 991-1110 (USA and Canada only) or 732 626-6466, option 1
 • Email interventional@angio.com • Visit www.AngioEduPRO.com, our educational website for healthcare professionals