

# Techna-Cut™ Biopsy Needles

The Techna-Cut™ Biopsy Needle is designed to perform simple, reliable biopsies of soft tissue. The Trocar Style needle tip allows for easy direct puncture while minimizing tissue damage.

- Numerically ordered centimeter markings on the outer cannula facilitate precise depth placement
- A precise cutting edge on the outer cannula contributes to a complete, intact core specimen
- A unique stylet design, combined with a locking syringe, provides vacuum for specimen retention in the outer cannula



BIOPSY

## Techna-Cut™ Biopsy Needle (Box of 10)

CATALOG NO.	DESCRIPTION
TCBN1615X	16G x 15 cm
TCBN1810X	18G x 10 cm
TCBN1815X	18G x 15 cm
TCBN1907X	19G x 7 cm
TCBN1910X	19G x 10 cm
TCBN1915X	19G x 15 cm
TCBN2107X	21G x 7 cm
TCBN2110X	21G x 10 cm
TCBN2115X	21G x 15 cm