

Access Needles

Stainless steel puncture needle with stylet, PTFE coated puncture needle with stylet

- Stainless steel puncture needle extra thin wall for maximum internal diameter.
- PTFE coating provides easy visualization during ultrasound guided needle insertion.

Stainless Steel Puncture Needle with Stylet **Box of 10**

- 5009 22 020** 22G (0.7 mm) x 20 cm, accepts guidewire .021" (0.53 mm)
- 5009 22 025** 22G (0.7 mm) x 25 cm, accepts guidewire .021" (0.53 mm)
- 5009 21 015** 21G (0.8 mm) x 15 cm, accepts guidewire .018" (0.46 mm)
- 5009 21 020** 21G (0.8 mm) x 20 cm, accepts guidewire .018" (0.46 mm)
- 5009 20 015** 20G (0.9 mm) x 15 cm, accepts guidewire .021" (0.53 mm)
- 5009 20 020** 20G (0.9 mm) x 20 cm, accepts guidewire .021" (0.53 mm)
- 5009 18 015** 18G (1.2 mm) x 15 cm, accepts guidewire .038" (0.97 mm)
- 5009 18 020** 18G (1.2 mm) x 20 cm, accepts guidewire .038" (0.97 mm)
- 5009 18 033** 18G (1.2 mm) x 33 cm, accepts guidewire .038" (0.97 mm)

PTFE Coated Puncture Needle

Box of 10

- 5015 18 015** 18G (1.2 mm) x 15 cm, accepts guidewire .038" (0.97 mm)
- 5015 18 020** 18G (1.2 mm) x 20 cm, accepts guidewire .038" (0.97 mm)

Trocar Needle

Box of 10

- 5007 21 015** 21G (1.2 mm) x 15 cm, accepts guidewire .018" (0.46 mm)
- 5007 18 015** 18G (1.2 mm) x 15 cm, accepts guidewire .035" (0.89 mm)
- 5007 18 020** 18G (1.2 mm) x 20 cm, accepts guidewire .035" (0.89 mm)

5009



5015



5007



Introducer Sheaths/Needles

Introducer Sheath/Needles

- Transparent or radiopaque sheaths for guidewire insertion or fluid aspiration.

Introducer Sheath Needle

Box of 10

- ISN 1915** 19G (1.0 mm) x 15 cm with 5F sheath, accepts guidewire .038" (0.97 mm)
- ISN 1920** 19G (1.0 mm) x 20 cm with 5F sheath, accepts guidewire .038" (0.97 mm)

ISN 1915

