

Hawkins™ Blunt Needle Access System

Hawkins™ Blunt Needle Access System

- Provides a lightweight method of targeting, remaining in place while position is confirmed. This is especially valuable in CT guided procedures
- The adjustable needle stop localizes the outer cannula to the target, preventing forward migration
- The outer cannula has an echogenic tip, promoting accurate placement under ultrasound guidance
- Numerically ordered centimeter markings on the outer cannula facilitate precise depth placement
- Accommodates most guidewires, biopsy needles, and biopsy instruments



Hawkins™ Blunt Needle Access System (Box of 5)

CATALOG NO.	DESCRIPTION
216057	16G x 5 cm
216107	16G x 10 cm
216157	16G x 15 cm
218057	18G x 5 cm
218107	18G x 10 cm
218157	18G x 15 cm
218207	18G x 20 cm