

Dilators

Radiopaque polyethylene dilators feature smooth rounded tips with gradual tapering facilitating tissue penetration.

- Perfect tapering improves tissue penetration.
- Luer lock hub for contrast injection and guidewire exchange.
- Radiopaque for fluoroscopic visualization.

Dilators

Box of 10

6104 35 020	4F x 20 cm, accepts guidewire .035" (0.89 mm)
6105 35 020	5F x 20 cm, accepts guidewire .035" (0.89 mm)
6105 38 020	5F x 20 cm, accepts guidewire .038" (0.97 mm)
6106 35 020	6F x 20 cm, accepts guidewire .035" (0.89 mm)
6106 38 020	6F x 20 cm, accepts guidewire .038" (0.97 mm)
6107 35 020	7F x 20 cm, accepts guidewire .035" (0.89 mm)
6107 38 020	7F x 20 cm, accepts guidewire .038" (0.97 mm)
6108 35 020	8F x 20 cm, accepts guidewire .035" (0.89 mm)
6108 38 020	8F x 20 cm, accepts guidewire .038" (0.97 mm)
6109 38 025	9F x 25 cm, accepts guidewire .038" (0.97 mm)
6110 38 025	10F x 25 cm, accepts guidewire .038" (0.97 mm)
6112 38 025	12F x 25 cm, accepts guidewire .038" (0.97 mm)
6114 38 025	14F x 25 cm, accepts guidewire .038" (0.97 mm)
6116 47 025	16F x 25 cm, accepts guidewire .047" (1.19 mm)



61xx