

Screw Dilators

A screw dilator atraumatically separates tissue, allowing dilation up to 14 French by using only one dilator.

- Minimal patient discomfort.
- Blunt dissection.
- Minimises risk of guidewire dislocation.
- Tissue dilation up to 14F using only one dilator.
- Radiopaque for fluoroscopic visualization.



A wide range of French sizes

Screw Dilators

Box of 10

6307 38 025 7F x 25 cm, accepts guidewire .038" (0.97 mm)

6308 38 025 8F x 25 cm, accepts guidewire .038" (0.97 mm)

6309 38 025 9F x 25 cm, accepts guidewire .038" (0.97 mm)

6310 38 025 10F x 25 cm, accepts guidewire .038" (0.97 mm)

6312 38 025 12F x 25 cm, accepts guidewire .038" (0.97 mm)

6314 38 025 14F x 25 cm, accepts guidewire .038" (0.97 mm)

6316 47 025 16F x 25 cm, accepts guidewire .047" (1.19 mm)

6318 52 025 18F x 25 cm, accepts guidewire .052" (1.32 mm)

63xx

