

SKATER® Drainage Catheters

7552/53

The SKATER® Drainage Catheter is made of soft polyurethane for optimal kink resistance and perfect configuration memory. The catheter is coated with HYDROFEEL® hydrophilic coating for minimal friction during insertion.

- Soft polyurethane for optimal kink resistance and perfect configuration memory.
- HYDROFEEL® coating on distal part for easy insertion.
- Large oval drainage holes and large lumen are designed for maximal drainage capacity.
- Minimal tissue trauma and ulcerations reduce patient discomfort.
- With locking or non-locking pigtail.



SKATER® Drainage Catheter • with locking pigtail **Box of 5**

- 7551 06 025** 6F x 25 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7551 06 035** 6F x 35 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7551 07 025** 7F x 25 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7551 07 035** 7F x 35 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7551 08 025** 8F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 08 035** 8F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 10 025** 10F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 10 035** 10F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 12 025** 12F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 12 035** 12F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 14 025** 14F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 14 035** 14F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7551 16 025** 16F x 25 cm, recommended guidewire .038" (0.97 mm), metal and flexible stiffening cannulae

SKATER® Drainage Catheter • with non-locking pigtail **Box of 5**

- 7511 06 025** 6F x 25 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7511 06 035** 6F x 35 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7511 07 025** 7F x 25 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7511 07 035** 7F x 35 cm, accepts guidewire .035" (0.89 mm), metal stiffening cannula
- 7511 08 025** 8F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 08 035** 8F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 10 025** 10F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 10 035** 10F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 12 025** 12F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 12 035** 12F x 35 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 14 025** 14F x 25 cm, accepts guidewire .038" (0.97 mm), metal and flexible stiffening cannulae
- 7511 16 025** 16F x 25 cm, recommended guidewire .038" (0.97 mm), metal and flexible stiffening cannulae