


SKATER™ All-Purpose Drainage Catheter

SKATER™ All-Purpose Drainage Catheters are suitable for a variety of applications. They feature a large catheter lumen, large drainage holes and a large lumen hub designed to promote maximum flow and drainage. SKATER™ is formulated for patient comfort, kink resistance, easy insertion and resistance to encrustation.

- Choice Lock trocar stylet allows for single-handed abscess drainage placement
- Hydrophilic coating allows for easy catheter insertion
- Locking shoulder tip ensures smooth insertion while preventing accordion effect
- Low profile, easy-to-use locking system with snap-off tab ensures secure pigtail deployment while maximizing patient comfort
- Polyurethane catheter material provides kink resistance
- Large skived drain holes promote maximum drainage while reducing catheter occlusion
- Includes SKATER™ All-Purpose Drainage Catheter with metal and flexible* stiffening cannulae, Choice Lock™ trocar stylet, removable pigtail straightener and locking system with snap-off tab



7565

SKATER™ All-Purpose Drainage Catheter w/Locking Pigtail (Box of 5)

CATALOG NO.	DESCRIPTION
756506020*	6F x 20 cm, accepts .035" guidewire
756507020*	7F x 20 cm, accepts .035" guidewire
756508025	8F x 25 cm, accepts .038" guidewire
756510025	10F x 25 cm, accepts .038" guidewire
756512025	12F x 25 cm, accepts .038" guidewire
756514025	14F x 25 cm, accepts .038" guidewire
756516025	16F x 25 cm, accepts .038" guidewire

SKATER™ All-Purpose Drainage Catheter w/Non-Locking Pigtail (Box of 5)

CATALOG NO.	DESCRIPTION
756006020*	6F x 20 cm, accepts .035" guidewire
756007020*	7F x 20 cm, accepts .035" guidewire



*Flexible stiffening cannula not available on 6F and 7F.

For more information, please call the Angiotech Customer Care Center at 877 991-1110 (USA and Canada only) or 732 626-6466, option 1
 • Email interventional@angio.com • Visit www.AngioEduPRO.com, our educational website for healthcare professionals