



SKATER™ Biliary Drainage Catheters

SKATER™ Biliary Drainage Catheters are used to drain bile from the bile duct and are placed using the Seldinger technique.

- Radiopaque marker band provides clear view of most proximal drain hole for accurate placement in the biliary tract
- Soft polyurethane provides kink-resistance and configuration memory
- Hydrophilic coating allows for easy insertion
- Large oval drainage holes and large lumen designed for maximum drainage capacity
- Includes SKATER™ Biliary Drainage Catheter with metal and flexible stiffening cannulae, removable pigtail straightener and locking system with snap-off tab



SKATER™ Biliary Drainage Catheter w/Locking Pigtail (Box of 5)

| CATALOG NO. | DESCRIPTION |
|-------------|--------------------------------------|
| 755308040 | 8F x 40 cm, accepts .038" guidewire |
| 755310040 | 10F x 40 cm, accepts .038" guidewire |
| 755312040 | 12F x 40 cm, accepts .038" guidewire |

For more information, please call the Angiotech Customer Care Center at 877 991-1110 (USA and Canada only) or 732 626-6466, option 1
 • Email interventional@angiotech.com • Visit www.AngioEduPRO.com, our educational website for healthcare professionals