

Connecting Tubes – Low Pressure

Non-leaking connection to standard luer lock connections

- Transparent soft tube.
- Standard luer lock connections.
- Wide range of lengths.
- Maximum pressure for tubes 280 psi.
- Maximum pressure for three-way stopcock 50 psi.

Connecting Tubes – Low Pressure

Box of 25

- 4010 04 050** 50 cm three-way stopcock/female-male, max pressure 50 psi
- 4010 11 010** 10 cm male-male luer lock, max pressure 280 psi
- 4010 11 025** 25 cm male-male luer lock, max pressure 280 psi
- 4010 12 010** 10 cm female-male luer lock, max pressure 280 psi
- 4010 12 025** 25 cm female-male luer lock, max pressure 280 psi
- 4010 12 050** 50 cm female-male luer lock, max pressure 280 psi
- 4010 12 080** 80 cm female-male luer lock, max pressure 280 psi
- 4010 13 025** 25 cm one-way stopcock/female-male luer lock, max pressure 280 psi
- 4010 13 050** 50 cm one-way stopcock/female-male luer lock, max pressure 280 psi
- 4010 13 080** 80 cm one-way stopcock/female-male luer lock, max pressure 280 psi
- 4010 43 010** Female/one-way stopcock – male luer, max pressure 280 psi

4010 11



4010 12



4010 13

