

HSG Catheters

HSG catheters are designed for diagnostic procedures, such as radiological hysterosalpingography or hysterosonography. The HSG catheter has a soft, concentric, non-latex balloon, which prevents leakage of saline or contrast media.

- Open ended, round distal tip provides atraumatic insertion with uniform contrast distribution
- The unique, malleable introduction sheath easily positions the catheter through the cervix into the uterine cavity
- Color coded stopcocks are attached to the catheter lumens; the blue stopcock is for balloon inflation and the white is for contrast media injection



6600

HSG Catheter (Box of 10)

| CATALOG NO. | DESCRIPTION |
|-------------|-------------------------------------|
| 660005040 | 5F x 40 cm, accepts .035" guidewire |
| 660007040 | 7F x 40 cm, accepts .038" guidewire |

HSG Procedure Trays (Box of 5)

| CATALOG NO. | DESCRIPTION |
|-------------|-------------|
| 660005000 | 5F x 40 cm |
| 660007000 | 7F x 40 cm |

Tray Contents Include:

- (1) HSG Catheter, 5F or 7F
- (1) 18" x 26" Fenestrated Drape
- (1) Welch Allyn® Speculum (compatible with light source #78810)
- (1) Intrauterine Sound
- (2) 3" x 3" Gauze Sponges
- (2) 8" Cotton Swabsticks
- (2) Wide Sponge Swabsticks
- (1) Inner Tray (may be used for Povidone Iodine Solution)
- (1) Povidone Iodine Ointment Packet – 1gm
- (1) Lubricating Jelly Packet – 5gm
- (2) Povidone Iodine Solution Packets – 1oz
- (1) Extension Tube Set – 30"
- (1) 20cc Luer Lock Syringe
- (1) 18G x 1– 1/2" Needle
- (1) Sanitary Napkin
- (1) 30" x 30" CSR Wrap