

Prostate Stabilization Needles

Prostate Stabilization Needles are designed to prevent movement of the prostate gland while either applicator-compatible or pre-loaded seeding needles are being inserted.

- Two-part hub design separates, locking the side barb in position
- The narrow side of the white, tear drop shaped hub indicates the deployment position of the side barb and the orientation of the needle point
- Each needle has an echogenic tip, promoting accurate placement under ultrasound

Prostate Stabilization Needles (Box of 20)

CATALOG NO.	DESCRIPTION
500118200	18G x 20 cm



Prostate Seeding Sets

Prostate Seeding Sets are designed to facilitate transperineal radioactive seed implant procedures.

- The 5003 series is designed for use with the Mick® 200-TP and 200-TPV applicators
- The 5005 needle cannula is more rigid for better trackability and is compatible with the Bard® EXPRESS Seeding Cartridge. It features one-half centimeter markings on the inner stylet to indicate the number of seeds in the needle after loading
- The 5006 needle cannula is more flexible for better steerability and features one-half centimeter markings on the inner stylet to indicate the number of seeds in the needle after loading



Prostate Seeding Sets (Box of 25)

CATALOG NO.	DESCRIPTION
500318200SX	For use with Mick® 200 – TP & 200 – TPV applicators 18G x 20 cm
500518200X	18G Rigid ETW x 20 cm
500618200X	18G Flexible UTW x 20 cm