

# Enteroclysis Catheters

The Enteroclysis Catheter is made of soft radiopaque polyvinylchloride for single or double contrast examination of the small bowel.

- 13F soft radiopaque catheter with a non-latex balloon.
- A large internal diameter allowing good contrast flow and a small external size for reduced patient discomfort.
- Stiff curved PTFE coated guidewire offers steerability of the catheter tip.
- The balloon prevents duodeno-gastric reflux of contrast media during examination.

## Enteroclysis Catheter

**Box of 2**

**6067 13 000** Set, 13F x 140 cm enteroclysis catheter with balloon, .052" (1.32 mm) x 150 cm stiffening wire with torque device

