

Duodenography Catheters and Sets

Radiopaque polyvinylchloride catheter with closed end and 6 side holes for even distribution of contrast. The catheter is fitted with female luer lock for secure connection to syringe. Available as a set or as single-packed products.

- Soft radiopaque flexible catheter for quick introduction.
- Round tip for atraumatic introduction.
- Malleable PTFE coated stiffening wire with knob to facilitate insertion and manipulation.



Duodenography Catheters

Box of 10

- 6006 10 150** 10F x 150 cm duodenography catheter
- 6006 12 125** 12F x 125 cm duodenography catheter
- 6006 12 150** 12F x 150 cm duodenography catheter
- 6006 14 125** 14F x 125 cm duodenography catheter
- 6006 14 150** 14F x 150 cm duodenography catheter

WORKER® Guidewires

Stiffening wire for duodenography

Box of 10

- 1175 52 150** .052" (1.32 mm) x 150 cm

Duodenography Sets

Box of 5

incl. Duodenography catheter and WORKER® Guidewire

- 6006 10 000** Set, 10F x 150 cm duodenography catheter, .052" (1.32 mm) x 150 cm stiffening wire with knob
- 6006 12 000** Set, 12F x 150 cm duodenography catheter, .052" (1.32 mm) x 150 cm stiffening wire with knob
- 6006 14 000** Set, 14F x 150 cm duodenography catheter, .052" (1.32 mm) x 150 cm stiffening wire with knob