

SKATER® Nephrostomy Kit with locking pigtail

- 18G puncture needle with stylet provides quick access for the strong, flexible and smooth WORKER® Amplatz wire.
- SKATER® Nephrostomy catheter with locking pigtail, and with metal or flexible stiffening cannula. HYDROFEEL® coating at the distal tip for easy insertion.
- Large lumen and large oval drainage holes designed for maximal drainage capacity.

SKATER® Nephrostomy Kit with non-locking pigtail


- 18G puncture needle with stylet provides quick access for the strong, flexible and smooth WORKER® Amplatz wire.
- SKATER® Nephrostomy catheter with non-locking pigtail, and with metal or flexible stiffening cannula. HYDROFEEL® coating at the distal tip for easy insertion.
- Large lumen and large oval drainage holes designed for maximal drainage capacity.

SKATER® Introducer Nephrostomy Kit


- Fine needle puncture with 21G trocar needle, .018" Access guidewire of either stainless steel or nitinol, and 6 French co-axial dilator provides atraumatic placement of .035" Amplatz guidewire.
- SKATER® Nephrostomy catheter with locking pigtail, and with metal or flexible stiffening cannula.
- Large oval drainage holes and large lumen are designed for maximal drainage capacity.
- Soft polyurethane for optimal kink resistance and perfect configuration memory.

SKATER® Introducer Biliary Drainage Kit with locking pigtail

- Fine needle puncture with 21G trocar needle, .018" Access guidewire of either stainless steel or nitinol, and 6 French co-axial dilator provides atraumatic placement of .035" Amplatz guidewire.
- SKATER® Biliary catheter with locking pigtail, and with metal or flexible stiffening cannula.
- Large oval drainage holes and large lumen are designed for maximal drainage capacity.
- Soft polyurethane for optimal kink resistance and perfect configuration memory.



Individual parts of the SKATER® Nephrostomy Kit ▲



SKATER® Introducer Biliary Drainage Kit with locking pigtail ▼

Box of 5

SKATER® Nephrostomy Kit • with locking pigtail

- 7573 06 600** 6F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F dilator, 6F x 35 cm SKATER® Nephrostomy catheter
- 7573 06 700** 6F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F dilator, 6F x 35 cm SKATER® Nephrostomy catheter, male/funnel connection for drainage bag
- 7573 07 000** 7F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F and 7F dilators, 7F x 35 cm SKATER® Nephrostomy catheter
- 7573 08 000** 8F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F and 8F dilators, 8F x 35 cm SKATER® Nephrostomy catheter
- 7573 08 600** 8F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F and 8F dilators, 8F x 35 cm SKATER® Nephrostomy catheter, male/funnel connection for drainage bag
- 7573 08 700** 8F Set, 22G (0.7 mm) x 20 cm chiba needle, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F and 8F dilators, 8F x 35 cm SKATER® Nephrostomy catheter, male/funnel connection for drainage bag
- 7573 10 000** 10F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F, 8F, and 10F dilators, 10Fx 35 cm SKATER® Nephrostomy catheter
- 7573 12 000** 12F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F, 8F, 10F, and 12F dilators, 12Fx 35 cm SKATER® Nephrostomy catheter
- 7573 14 700** 14 Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F, 8F, 10F, 12F, and 14F dilators, 14Fx 35 cm SKATER® Nephrostomy catheter

Box of 5

SKATER® Nephrostomy Kit • with non-locking pigtail

- 7503 08 000** 8F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F and 8F dilators, 8Fx 35 cm SKATER® Nephrostomy catheter
- 7503 08 100** 8F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 100 cm 3 mm J HD guidewire, 6F and 8F dilators, 8Fx 35 cm SKATER® Nephrostomy catheter
- 7503 10 000** 10F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F, 8F, and 10F dilators, 10Fx 35 cm SKATER® Nephrostomy catheter
- 7503 12 000** 12F Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F, 8F, 10F, and 12F dilators, 12Fx 35 cm SKATER® Nephrostomy catheter
- 7503 14 000** 14 Set, 18G (1.2 mm) x 15 cm puncture needle with stylet, .035" (0.89 mm) x 90 cm Amplatz guidewire, 6F, 8F, 10F, 12F, and 14F dilators, 14Fx 35 cm SKATER® Nephrostomy catheter

Box of 5

SKATER® Introducer Nephrostomy Kit • with locking pigtail

- 7673 08 000** 8F Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm stainless steel access wire, 6F x 18 cm co-axial dilator, .035" (0.89 mm) x 90 cm Amplatz guidewire, 8Fx 20 cm dilator, 8F x 35 cm catheter with locking pigtail
- 7673 08 300** 8F Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access wire, 6F x 18 cm co-axial dilator, .035" (0.89 mm) x 90 cm Amplatz guidewire, 8Fx 20 cm dilator, 8F x 35 cm catheter with locking pigtail
- 7673 10 000** 10F Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm stainless steel access wire, 6F x 18 cm co-axial dilator, .035" (0.89 mm) x 90 cm Amplatz guidewire, 8F and 10F dilators, 10F x 35 cm catheter with locking pigtail
- 7673 10 300** 10F Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access wire, 6F x 18 cm co-axial dilator, .035" (0.89 mm) x 90 cm Amplatz guidewire, 8F and 10F dilators, 10F x 35 cm catheter with locking pigtail
- 7673 12 000** 12F Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm stainless steel access wire, 6F x 18 cm co-axial dilator, .035" (0.89 mm) x 90 cm Amplatz guidewire, 8F, 10F, and 12F dilators, 12F x 35 cm catheter with locking pigtail

Box of 5

SKATER® Introducer Biliary Drainage Kit • with locking pigtail

- 7653 08 300** 8F Set, 21G (0.8 mm) x 15 cm trocar needle, .018" (0.46 mm) x 80 cm nitinol access wire, 6F x 26 cm co-axial dilator, .035" (0.89 mm) x 150 cm Amplatz guidewire, 8F x 20 cm dilator, 8F x 40 cm catheter with locking pigtail
- 7653 10 300** 10F Set, 21G x 15 cm trocar needle, .018" x 80 cm nitinol access wire, 6F x 26 cm co-axial dilator, .035" (0.89 mm) x 150 cm Amplatz guidewire, 8F x 20 cm and 10Fx 20 cm dilators, 10F x 40 cm catheter with locking pigtail