

# Multiband Ligator

## Barrel with 6 bands

**Used for endoscopic ligation of Esophageal varices at or above gastroesophageal junction**

Fits on wide range of endoscopes  
Second last Band light colour

**Standard Band Ligator**  
requires 9.5 to 11.2 mm  
endoscope diameter

**Paediatric Band Ligator**  
requires 8.0 to 9.5 mm  
endoscope diameter



|            | Art. No. | Product name                                                   |
|------------|----------|----------------------------------------------------------------|
| Standard   | SGK-6    | Set incl. Applicator barrel, thread puller and flushing needle |
|            | SGB-6    | Barrel and thread puller only                                  |
| Paediatric | SGK-6-S  | Set incl. Applicator barrel, thread puller and flushing needle |
|            | SGB-6-S  | Barrel and thread puller only                                  |