


Shark Tooth Forcep Alligator Forcep

Used to retrieve foreign body from Gastrointestinal Tract

Shark Tooth Forcep can be used with ERCP scope to retrieve Biliary-Stent



Shark Tooth Forcep

Alligator Forcep

Art. No.	Product name	Outer Diameter	Working length	Working Channel Ø
SSTF-230	Shark Tooth Forcep Rat tooth jaws	2.5 mm	230 cm	2.8 mm
SATF-170	Alligator Forcep Alligator jaws	2.5 mm	170 cm	2.8 mm