


Retrieval Spiral Basket

Nitinol

Used to retrieve foreign bodies from Gastrointestinal Tract as well as Pancreatic Stones with Spiral Basket through the scope

Not for Lithotripsy



Art. No.	Basket Ø	Working Length	Sheath Ø	Working Channel Ø
WS-1819GL520	20 mm	190 cm	1.8 mm	2.0 mm
WS-2419GL520	20 mm	190 cm	2.4 mm	2.8 mm
WS-2415GL530	30 mm	150 cm	2.4 mm	2.8 mm