

# Monofilament Memory Baskets

## Wire Guided

**Used to retrieve foreign bodies from Gastrointestinal Tract through the scope**

Handle light weight with flush port for Contrast-Drug-Flush  
Retains shape for longer time  
Centrally located wireguide entrance

Working length 230 cm  
Compatible with working Channel 2.8 mm diameter

**Not for Lithotripsy**



Art. No.	Product name	Working Channel Ø	Basket Ø
SMMB-WG-30	Monofilament Memory Basket Wire Guided	2.8 mm	30 mm
SMMB-WG-35	Monofilament Memory Basket Wire Guided	2.8 mm	35 mm
SMMB-WG-40	Monofilament Memory Basket Wire Guided	2.8 mm	40 mm