


Self expandable Esophageal Stent

Used to treat Tracheoesophageal Fistulas and expansion treatment of the Esophageal

Stenosis caused by Esophagael Carcinoma, Chemical Injuries, Cardiac Carcinoma, Anastomatic stenosis, and other Traumas, which are surgical contraindicated.

Made of NITI – Nickel/Titanium – radiopac – Stent silicone coated with 18 mm diameter Introducer System



Art. No.	Product name	Stent Length	Stent
S-ES-5	Niti – Stent	50 mm long	Silicone coated
S-ES-8	Niti – Stent	80 mm long	Silicone coated
S-ES-10	Niti – Stent	100 mm long	Silicone coated
S-ES-12	Niti – Stent	120 mm long	Silicone coated
S-ES-15	Niti – Stent	150 mm long	Silicone coated