

# Esophageal Balloon Dilator

**Used to dilate strictures of the Gastrointestinal Tract including strictures of the Esophagus through the scope**

Shaft Ø 2.5 mm

Guidewire compatibility 0.035"



Art. No.	Product Name	Balloon length	Balloon Ø	Working length
EBD-0630	Esophageal Balloon Dilator	30 mm	6 mm	180 cm
EBD-0830	Esophageal Balloon Dilator	30 mm	8 mm	180 cm
EBD-1040	Esophageal Balloon Dilator	40 mm	10 mm	180 cm
EBD-1280	Esophageal Balloon Dilator	80 mm	12 mm	180 cm
EBD-1480	Esophageal Balloon Dilator	80 mm	14 mm	180 cm
EBD-1680	Esophageal Balloon Dilator	80 mm	16 mm	180 cm
EBD-1880	Esophageal Balloon Dilator	80 mm	18 mm	180 cm
EBD-18230	Esophageal Balloon Dilator	80 mm	18 mm	230 cm
EBD-20230	Esophageal Balloon Dilator	80 mm	20 mm	230 cm

other sizes on request