

Achalasia-Cardia Balloon Dilator

**Used for dilation of the lower esophageal sphincter
in patient with achalasia (stomach entry)**

High Radial Vector Force
Atraumatic Tip
Radiopaque marker inside the balloon

Shaft Ø 5 mm
Guidewire compatibility 0.038"



Art. No.	Product name	Balloon Ø	Balloon length	Working length
SABD-25	Achalasia Balloon Dilator	25 mm	80 mm	110 cm
SABD-30	Achalasia Balloon Dilator	30 mm	80 mm	110 cm
SABD-35	Achalasia Balloon Dilator	35 mm	80 mm	110 cm